

# Automatic Densimeters



## INTRODUCTION

The Automatic Densimeter offers high precision along with ease of operation. An electronic balance determines specific gravity with the single push of a button. An optional attachment allows measurement of the specific gravity of liquids or the rate of change of volume for rubber and other materials

## APPLICATIONS

**Rubber, Plastic, Metals, Glass, Ceramic, Food,  
Wood, Medical**

## SPECIFICATIONS

	Model D-H 300	Model D-H 100
Resolution	.001g	.0001g
Capacity	300 g	100 g

## FEATURES

- Measures specific gravity of a wide range of substances
- Specific gravity range 0.1 to 2.0
- Only 20 seconds to measure one sample
- Simple operation —mount sample, push start button, dismount sample
- Distilled water is used as standard liquid, but other liquids may be used
- Many attachments available for versatility in testing

## INSTRUMENT SIZE

Depth:	420 mm (17 in.)
Height:	890 mm (35 in.)
Width:	460 mm (18 in.)

## WEIGHT

30 kg (66 lb)

## CONNECTIONS

Electrical: Specify voltage requirements when ordering

## ORDERING INFORMATION

Model No. D-H 300, catalog number 21-09-01

Model No. D-H 100, catalog number 21-09-02

**The TMI Group  
of Companies**



[www.testingmachines.com](http://www.testingmachines.com)

**Testing Machines Inc.** 40 McCullough Drive, New Castle, DE 19720 USA

Tel: 800-678-3221 (USA) or (302) 613-5600, Fax: (302) 613-5619, e-mail: [info@testingmachines.com](mailto:info@testingmachines.com)

**Büchel BV t/a Messmer Büchel** Fokkerstraat 24, 3905 KV Veenendaal, Netherlands

Tel: +31 (0)318 521500, Fax: +31 (0)318 540358, e-mail: [buchel@testingmachines.com](mailto:buchel@testingmachines.com)

**Lako Tool** 7400 Ponderosa Road, Perrysburg, OH 43552 USA

Tel: (419) 662-5256, Fax: (419) 662-8225

**Lawson Hemphill** 1658 G A R Highway, Swansea, MA 02777 USA

Tel: (508) 679-5364, Fax: (508) 679-5396, e-mail: [information@lawsonhemphill.com](mailto:information@lawsonhemphill.com)

**Adamel Lhomargy SAS Z.A.** de l'Habitat, Batiment 6, Route d'Ozoir, 77680, Roissy en Brie, France

Tel: +33 (0)1 64409210, Fax: +33 (0)1 64409211, e-mail: [adamellhomargy@testingmachines.com](mailto:adamellhomargy@testingmachines.com)

**TMI Canada** 40 McCullough Drive, New Castle, DE 19720 USA

Tel: (302) 613-5600, Fax: (302) 613-5619 e-mail: [canada@testingmachines.com](mailto:canada@testingmachines.com)