

Automatic Densimeters



SPECIFICATIONS

	Model D-H 300	Model D-H 100
Resolution	.001g	.0001g
Capacity	300 g	100 g

FEATURES

- Measures specific gravity of a wide range of substances
- Specific gravity range 0.1 to 2.0
- Only 20 seconds to measure one sample
- Simple operation —mount sample, push start button, dismount sample
- Distilled water is used as standard liquid, but other liquids may be used
- Many attachments available for versatility in testing

INSTRUMENT SIZE

Depth:	420 mm (17 in.)
Height:	890 mm (35 in.)
Width:	460 mm (18 in.)

WEIGHT

30 kg (66 lb)

CONNECTIONS

Electrical: Specify voltage requirements when ordering

ORDERING INFORMATION

Model No. D-H 300, catalog number 21-09-01

Model No. D-H 100, catalog number 21-09-02

INTRODUCTION

The Automatic Densimeter offers high precision along with ease of operation. An electronic balance determines specific gravity with the single push of a button. An optional attachment allows measurement of the specific gravity of liquids or the rate of change of volume for rubber and other materials

APPLICATIONS

Rubber, Plastic, Metals, Glass, Ceramic, Food, Wood, Medical

The TMI Group of Companies



www.testingmachines.com

Testing Machines Inc. 40 McCullough Drive, New Castle, DE 19720 USA

Tel: 800-678-3221 (USA) or (302) 613-5600, Fax: (302) 613-5619, e-mail: info@testingmachines.com

Messmer Instruments Ltd. Unit F1 Imperial Business Estate, West Mill, Gravesend, Kent DA11 0DL UK

Tel: +44 (0)1474 566488, Fax: +44 (0)1474 560310, e-mail: messmer@testingmachines.com

Büchel BV Fokkerstraat 24, 3905 KV Veenendaal, Netherlands

Tel: +31 (0)318 521500, Fax: +31 (0)318 540358, e-mail: buchel@testingmachines.com

Adamel Lhomargy SAS Z.A. de l'Habitat, Batiment 6, Route d'Ozoir, 77680, Roissy en Brie, France

Tel: +33 (0)1 64409210, Fax: +33 (0)1 64409211, e-mail: adamellhomargy@testingmachines.com

TMI Canada 6 Ronald Drive, Montreal West, Quebec H4X 1M8

Tel: (514) 489-8491, Fax (514) 489-6044, e-mail: canada@testingmachines.com