

Valley Laboratory Beater



INTRODUCTION

The Valley Laboratory Beater is used to beat pulp in a uniform, standard and reproducible way. A single 360 gram sample of pulp can be beaten in the large (25 liter) tub. Standard procedure for grinding and calibrating will ensure reproducible results.

APPLICATIONS

Pulp

SPECIFICATIONS

- 25 liter tub
- Conforms to National and International Standards
ISO 5264/1, TAPPI T200, BS 6094 Part 1, SCAN C25

FEATURES

- Large pulp sample allows for progressive investigation of work done from one single charge
- Adjustable and repeatable beating pressure is produced by a weighted lever
- Sturdy construction

INSTRUMENT SIZE

Depth: 650 mm (25.5 in.)
Height: 970 mm (38 in.)
Width: 950 mm (37.5 in.)

WEIGHT

approx. 400 kg (875 lb)

CONNECTIONS

Electrical: Specify voltage requirements when ordering

ORDERING INFORMATION

Catalog Number	Material
72-13-01	Cast Iron
72-13-02	Cast Bronze
72-16	Stainless Steel

The TMI Group of Companies



www.testingmachines.com

Testing Machines Inc. 2 Fleetwood Court, Ronkonkoma, NY 11779 USA
Tel: 800-678-3221 (USA) or (631) 439-5400, Fax: (631) 439-5420, e-mail: info@testingmachines.com

Messmer Instruments Ltd. Unit F1 Imperial Business Estate, West Mill, Gravesend, Kent DA11 0DL UK
Tel: +44 (0)1474 566488, Fax: +44 (0)1474 560310, e-mail: messmer@testingmachines.com

Büchel BV Fokkerstraat 24, 3905 KV Veenendaal, Netherlands
Tel: +31 (0)318 521500, Fax: +31 (0)318 540358, e-mail: buchel@testingmachines.com

Adamel Lhomargy SAS Z.A. de l'Habitat, Batiment 6, Route d'Ozoir, 77680, Roissy en Brie, France
Tel: +33 (0)1 64409210, Fax: +33 (0)1 64409211, e-mail: adamellhomargy@testingmachines.com

TMI Canada 6 Ronald Drive, Montreal West, Quebec H4X 1M8
Tel: (514) 489-8491, Fax (514) 489-6044, e-mail: canada@testingmachines.com