

Pulp Disintegrator



INTRODUCTION

Disintegration of a pulp sample is the mechanical treatment, in water, of interlaced fibers to separate them from one another.

APPLICATIONS

Pulp

SPECIFICATIONS

- Meets TAPPI T-205, CPPA C.4 and ISO 5263

FEATURES

- Digital counter
- Includes lightweight (1.1 kg) polymer container
- 316 Stainless cover
- Integrated sensors for container proximity and cover lock enhance operator safety

OPTIONS

- Hot disintegrator-heated stock container maintains temperature of pre-heated pulp samples

INSTRUMENT SIZE

Depth: 521 mm (20.5 in.)
Height: 425 mm (16.75 in.)
Width: 305 mm (12 in.)

WEIGHT

40 kg (87 lb)

CONNECTIONS

Electrical: Specify voltage requirements when ordering

ORDERING INFORMATION

Pulp Disintegrator	Catalog number 73-18
Hot Disintegrator (w/o container)	Catalog number 73-19
3 liter heated container	Catalog number 73-19-01
7 liter container	Catalog number 73-19-02

The TMI Group of Companies



www.testingmachines.com

Testing Machines Inc. 2 Fleetwood Court, Ronkonkoma, NY 11779 USA
Tel: 800-678-3221 (USA) or (631) 439-5400, Fax: (631) 439-5420, e-mail: info@testingmachines.com

Messmer Instruments Ltd. Unit F1 Imperial Business Estate, West Mill, Gravesend, Kent DA11 0DL UK
Tel: +44 (0)1474 566488, Fax: +44 (0)1474 560310, e-mail: messmer@testingmachines.com

Büchel BV Fokkerstraat 24, 3905 KV Veenendaal, Netherlands
Tel: +31 (0)318 521500, Fax: +31 (0)318 540358, e-mail: buchel@testingmachines.com

Adamel Lhomargy SAS Z.A. de l'Habitat, Batiment 6, Route d'Ozoir, 77680, Roissy en Brie, France
Tel: +33 (0)1 64409210, Fax: +33 (0)1 64409211, e-mail: adamellhomargy@testingmachines.com

TMI Canada 6 Ronald Drive, Montreal West, Quebec H4X 1M8
Tel: (514) 489-8491, Fax (514) 489-6044, e-mail: canada@testingmachines.com